



MOMENTUM

EQUINE

SOLUTIONS

Revolutionary Imaging Systems




[2020imaging.net/equine](http://2020imaging.net/equine)





# Momentum EQ All In One Portable Solution





 Lightweight Design Under 18 lbs


 Counterbalance for stability


 Ergonomic Panel and Tablet Holders

 Allows for easy handling


 Custom Tablet Enclosure



 All in one charging


 Requires one power cable



 Rugged Lightweight Aluminum Design

 Scratch Resistant Finish

# 10x 12 Premium DR Panel



Lightest panel & holder at 4lbs	16-bit AD Conversion
Scintillator: Csl (Cesium)	Battery Life: > 150 images/4hrs
Cassette sized	Rapid Charge enabled
Imaging Area: 10" x 12"	Warranty: 3 years
Integrated Cable Management	7 year lithium-ion capacitor

# Veterinary Specific Software



- |                                |                           |                   |                               |
|--------------------------------|---------------------------|-------------------|-------------------------------|
| • All-in-One Software Solution | • User-friendly Interface | • Clinical Review | • Automated Image Processing  |
| • Equine Specific Protocols    | • DR Console              | • Image Storage   | • Produces an optimized image |

Advanced · Superior · Efficient · Reliable · Affordable



[2020imaging.net/equine](http://2020imaging.net/equine)