



DIGITAL X-RAY

VETERINARY

SOLUTIONS

Revolutionary Imaging Systems




2020imaging.net/vet



Wireless, Fixed, & Hybrid DR



MOMENTUM 14X17 WIRELESS FLATPANEL



- Wireless Interface
- Rapid Charge
- Cesium Scintillator
- Lithium Ion Capacitor
- Featherlight 5.5 lbs
- Cassette Sized
- 16-bit Images

HYBRID 17X17 FLATPANEL



- Wireless or Tethered Interface
- Cesium Scintillator
- Lithium Ion Battery
- Lightweight 9.5 lbs
- Cassette Sized
- 16-bit Images

WIRELESS 14X17 FLATPANEL

- Wireless Interface
- Cesium Scintillator
- Lithium Ion Battery
- Lightweight 8 lbs
- Cassette Sized
- 16-bit Images



TETHERED 17X17 DETECTOR **NEW**

- Tethered Interface
- Cesium Scintillator
- 16-bit Images
- Cassette Sized



Retrofit and New X-Ray



* RETROFIT ANY X-RAY

Wireless & Fixed DR Compatible



VXR X-RAY SYSTEM

- Superior Image Quality
- Anatomical Programmed
- Table-side Acquisition
- Acquisition Touch-Screen
- Fixed 40" SID integrated tube
- 32" x 60" 4-way float table top
- Exposure Foot Switch Pedal
- 32kW, 500mAs, 40-125kVp

2020imaging.net/straight-arm

Veterinary Specific Software

ACQUISITION INTERFACE



- Under 6 second ⌚ Processing
- Background Submission
- Preset X-ray Exams
- Animal Specific
- Anatomical Program (APR)
- Fully Customizable
- Advanced Image Processing
- Image Adjustment Tools
- Touchscreen Enabled
- User Friendly Interface

☰ Multi-funtion Studylist ⬇️ 🗨️ Advanced Viewer ⬆️

- PM Software Integration**
- Burn/Import CDs
- Custom Search/Query
- Customized User Settings
- Vertibral Heart Score Tool incl.
- Fast, User Friendly Interface
- Image Fine Tuning
- Exam Room Viewing



MOBILE VIEWING WITH M-OPAL[†]

iPad®, iPhone®, Android™, Blackberry®, any mobile or smart device!



- Lightweight & Fast
- Patient List & Viewer
- Simple Annotations Enabled
- Tap, Drag, and Pinch
- Patient Filter & Search



CLOUD-BACKUP INCLUDED^{††}

- Live Synchronization
- Fully Monitored
- 448-bit Encryption (traffic)
- Fully Customizable
- Enterprise Grade
- Recover Deleted Files



Advanced
Superior
Efficient
Reliable
Affordable

2020imaging.net/vet

 **20/20 Imaging**
a division of Konica Minolta Healthcare