

ELECTROTHERAPY | ULTRASOUND | COMBINATION

25

S E R I E S




DYNATRON® 25 SERIES™

3 & 5 CHANNEL MODELS with UNEQUALED COMBOPLUS™ COMBINATION THERAPY

The only 3 and 5 channel devices available. All channels allow fully-independent treatment setups offering Interferential, Premodulated, High Volt, Biphasic, Russian and Fixed Frequency IFC/Premod. In addition, the 925 and 825 include Dynatronics' Ultrasound Comboplus feature with the power to deliver up to 5 channels of stim and Ultrasound — all at the same time.





3 FREQUENCY ULTRASOUND

Dynatron 25 Series models offer 1, 2, and 3 MHz frequencies for the greatest flexibility in depth of treatment. Choose 1 MHz for deep treatments, 2 MHz for moderate depth, or 3 MHz for superficial depth.



TARGET & TARGET SWEEP for INTERFERENTIAL TREATMENTS

There's no easier, more efficient way to focus treatment precisely where it's needed! Simply glide your finger across the TARGET touchpad to move the center of interference to the site of your patient's pain. Available on all 25 Series electrotherapy devices.



IMPRESSIVELY SIMPLE

Even with nearly unlimited flexibility, setting up a treatment has never been easier. Most treatments can be set up in as few as three key presses. Treatment parameters are highly visible, easy to read, and changed at the touch of a fingertip.



Dynatron 925
5 Channel Combo Unit



Dynatron 825
3 Channel Combo Unit



Dynatron 125
Ultrasound Unit



Dynatron 625
5 Channel Stim Unit



Dynatron 525
3 Channel Stim Unit

	925	825	625	525	125	
ELECTROTHERAPY	IFC	✓	✓	✓	✓	—
	Premod	✓	✓	✓	✓	—
	Fixed Frequency IFC/Premod	✓	✓	✓	✓	—
	Biphasic	✓	✓	✓	✓	—
	Russian	✓	✓	✓	✓	—
	High Volt	✓	✓	✓	✓	—
	Combo Electrotherapy / Ultrasound	✓	✓	—	—	—
	Ultrasound	✓	✓	—	—	✓
	Available Electrotherapy Channels	5	3	5	3	0



7030 Park Centre Drive | Salt Lake City, UT | 84121

toll-free [800] 874-6251 | local [801] 568-7000 | fax [801] 568-7711

www.dynatronics.com | email: info@dynatron.com