

AMERICOMP AC1 CHIROPRACTIC SYSTEM

INCLUDING SPECTRA HF30 HIGH FREQUENCY GENERATOR

PRICING STARTS AT \$17,800.00

SPECTRA HF30 HIGH FREQUENCY GENERATOR WITH THE FOLLOWING:

A3306 Tubestand

- Floor-to-wall tubestand with 5 ft. rail
- Mechanical locks
- Angulation dial and operator handgrips
- Angulating platform tube mount with 10" arm

A6000 Wallstand Package

- Counterbalanced column, floor-to-wall mounted
- Grid cabinet, 17" x 17"
- Grid: 8:1 ratio, 103 LPI
- Regular-duty cassette tray
- Mechanical locks

X-Ray Tube

- 1.0 – 2.0mm focal spots, 140,000 heat units with 180 degree horns

C600-15 Cables

- 15 ft. high voltage cables with federal terminations

D800-1 Collimator

- Certified manual collimator with light field, non-swivel