

AMRAD MEDICAL AE4 STATIONARY X-RAY SYSTEM

INCLUDING HIGH FREQUENCY GENERATOR

PRICING STARTS AT \$23,300.00

AMRAD® Medical AE4 303E/301E Generator with following:

S200 Stationary Table

- Stationary-top radiographic table, 76" x 30"
- Welded base construction
- Enclosed pedestal
- Grid cabinet with 50" travel
- Grid: 10:1 ratio, 103 LPI
- Heavy-duty cassette tray 03092
- 300 lb. rated patient load

S105 Tubestand

- Floor-to-wall or floor to ceiling tubestand with 8 ft. rail
- Electric locks
- Angulation dial and operator handgrips
- Angulating tube mount

S109 Wallstand Package

- Counterbalanced column, floor-to-wall mounted
- Grid Cabinet, 17 in. x 17 in.
- Grid: 10:1 ratio, 103 LPI
- Heavy-duty cassette tray 03092
- Specify left or right load
- Electric locks

Standard X-Ray Tube (Specified by kW rating)

- 1.0 – 2.0 mm focal spots, 140,000 heat units (30 or 32 kW generator) – C458

Cables J995

- 25 ft. High Voltage cables, 150 kVp

07385 Collimator Ralco 221

- Certified 150 kVp manual collimator with swivel mount
- LED light field and tape measure
- Transverse laser alignment to receptor tray