DEL MEDICAL

RT100

FOUR-WAY FLOAT TOP TABLE

This sturdy table is ideally suited for a variety of general radiographic needs. Its compact size make the RT100 easy to use in small hospitals, clinics, ambulatory care centers and private offices. The four-way floating table top glides freely on a precision roller bearing system and can support patients up to 700 lbs. (318 kg)

Key Highlights

Four-way Float Top

- Glides effortlessly in all directions

Heavy-duty Electric Table Locks with Foot Treadle Release

700 lbs. Patient Weight Capacity

All Steel Structure Welded Base

- Provides maximum strength and reliability

Pedestal Base Can House Generator Electronics to Save Space

Smooth, Stain Resistant Fiber Resin Top



 $\ensuremath{^*}$ All specifications subject to change without notice

RT100 Four-Way Float Top Table	
Load capacity:	700 lbs. (318 kg)
Table base length:	43" (109 cm)
Table base width:	29" (74 cm)
Tabletop movement:	Longitudinal, 51.50" (131 cm) Transverse, 9.5" (24 cm)
Tabletop-to-film distance:	3.25" (8 cm)
Total bucky travel:	Longitudinal 11.25" (28 cm)
Tabletop length:	86.66" (220 cm)
Tabletop width:	32" (81 cm)
Weight:	387 lbs. (176 kg)
Power requirements:	120/230-240 VAC
Options:	78" (198 cm) tabletop length / Variety of grid options / Lateral cassette holder / Table-side patient handgrips / Compression band





