

AMRAD MEDICAL ELEVATING AE1 X-RAY SYSTEM SYSTEM

INCLUDING HIGH FREQUENCY GENERATOR—

PRICING STARTS AT \$29,199.00

AMRAD® Medical AE1 303E/301E Generator with following:

S211 Elevator / Four-Way Float-Top Table

- 84 in. x 30 in. flat, seamless table top
- Variable table height, 23 in. to 35 in.
- +/- 5 in. transverse travel
- 32 in. longitudinal travel
- Foot treadle operation of table top and height movements
- Grid Cabinet with electric lock
- Grid: 10:1 ratio, 103 LPI
- Heavy-duty cassette tray 03092
- Operates on 120VAC
- 400 lb. rated patient load

S225 Transverse Rotational Tubestand

- Floor-to-wall or floor-to-ceiling mount with 8 ft. rails
- Rotational base, +/- 360° with 90° detents
- Transverse tubearm +/- 5 in., platform tube mount
- Angulation dial, operator handgrips, and electric locks

S109 Wallstand Package

- Counterbalanced column, floor-to-wall mounted
- Grid Cabinet, 17 in. x 17 in.
- Grid: 10:1 ratio, 103 LPI
- Heavy-duty cassette tray 03092
- Specify left or right load
- Electric locks

Standard X-Ray Tube

- 1.0 – 2.0 mm focal spots, 140,000 heat units

Cables J995

- 25 ft. High Voltage cables, 150 kVp

07385 Collimator Ralco 221

- Certified 150 kVp manual collimator with swivel mount
- LED light field and tape measure
- Transverse laser alignment to receptor tray